


(6-Bromopyridin-3-yl)methanol

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 122306-01-8 |
| Formula: | C ₆ H ₆ BrNO |
| Molecular Weight: | 188.022 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|-------------|---|
| Description | (6-Bromopyridin-3-yl)methanol is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery. |
|-------------|---|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 5.3186 mL | 26.5929 mL | 53.1858 mL |
| 5 mM | 1.0637 mL | 5.3186 mL | 10.6372 mL |
| 10 mM | 0.5319 mL | 2.6593 mL | 5.3186 mL |
| 50 mM | 0.1064 mL | 0.5319 mL | 1.0637 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481