

## 1,2-Dianilinoethane


## Chemical Properties

CAS No. : 150-61-8

Formula: C<sub>14</sub>H<sub>16</sub>N<sub>2</sub>

Molecular Weight: 212.2902

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1,2-Dianilinoethane is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7105 mL	23.5527 mL	47.1054 mL
5 mM	0.9421 mL	4.7105 mL	9.4211 mL
10 mM	0.4711 mL	2.3553 mL	4.7105 mL
50 mM	0.0942 mL	0.4711 mL	0.9421 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481