


3-Acetyl-5-bromopyridine

Chemical Properties

CAS No. :	38940-62-4
Formula:	C ₇ H ₆ BrNO
Molecular Weight:	200.033
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	3-Acetyl-5-bromopyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9993 mL	24.9963 mL	49.9925 mL
5 mM	0.9999 mL	4.9993 mL	9.9985 mL
10 mM	0.4999 mL	2.4996 mL	4.9993 mL
50 mM	0.100 mL	0.4999 mL	0.9999 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481