


3-Nitrosalicylic acid

Chemical Properties

CAS No. :	85-38-1
Formula:	C7H5NO5
Molecular Weight:	183.1183
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	3-Nitrosalicylic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4609 mL	27.3045 mL	54.609 mL
5 mM	1.0922 mL	5.4609 mL	10.9218 mL
10 mM	0.5461 mL	2.7304 mL	5.4609 mL
50 mM	0.1092 mL	0.5461 mL	1.0922 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481