


## 2,5-Dichloropyridine

## Chemical Properties

CAS No. :	16110-09-1
Formula:	C <sub>5</sub> H <sub>3</sub> Cl <sub>2</sub> N
Molecular Weight:	147.99
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	2,5-Dichloropyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.7572 mL	33.7861 mL	67.5721 mL
5 mM	1.3514 mL	6.7572 mL	13.5144 mL
10 mM	0.6757 mL	3.3786 mL	6.7572 mL
50 mM	0.1351 mL	0.6757 mL	1.3514 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481