


## 2-Amino-3-chlorobenzoic acid

## Chemical Properties

CAS No. :	6388-47-2
Formula:	C <sub>7</sub> H <sub>6</sub> ClNO <sub>2</sub>
Molecular Weight:	171.581
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	2-Amino-3-chlorobenzoic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8282 mL	29.1409 mL	58.2819 mL
5 mM	1.1656 mL	5.8282 mL	11.6564 mL
10 mM	0.5828 mL	2.9141 mL	5.8282 mL
50 mM	0.1166 mL	0.5828 mL	1.1656 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481