


2-Bromo-4-chloropyridine

Chemical Properties

CAS No. :	22918-01-0
Formula:	C ₅ H ₃ BrClN
Molecular Weight:	192.441
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-Bromo-4-chloropyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1964 mL	25.9821 mL	51.9642 mL
5 mM	1.0393 mL	5.1964 mL	10.3928 mL
10 mM	0.5196 mL	2.5982 mL	5.1964 mL
50 mM	0.1039 mL	0.5196 mL	1.0393 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481