


## 4-Amino-2-chloropyrimidine

## Chemical Properties

CAS No. :	7461-50-9
Formula:	C <sub>4</sub> H <sub>4</sub> ClN <sub>2</sub>
Molecular Weight:	129.548
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	4-Amino-2-chloropyrimidine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.719 mL	38.5951 mL	77.1903 mL
5 mM	1.5438 mL	7.719 mL	15.4381 mL
10 mM	0.7719 mL	3.8595 mL	7.719 mL
50 mM	0.1544 mL	0.7719 mL	1.5438 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481