


2-Chloronicotinaldehyde

Chemical Properties

CAS No. :	36404-88-3
Formula:	C ₆ H ₄ ClNO
Molecular Weight:	141.555
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-Chloronicotinaldehyde is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.0641 mL	35.3207 mL	70.6414 mL
5 mM	1.4128 mL	7.0641 mL	14.1283 mL
10 mM	0.7064 mL	3.5321 mL	7.0641 mL
50 mM	0.1413 mL	0.7064 mL	1.4128 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481