


## 5-Fluoro-2-methoxyphenylboronic acid

## Chemical Properties

CAS No. :	179897-94-0
Formula:	C <sub>7</sub> H <sub>8</sub> BF <sub>3</sub> O <sub>3</sub>
Molecular Weight:	169.946
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-Fluoro-2-methoxyphenylboronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8841 mL	29.4204 mL	58.8408 mL
5 mM	1.1768 mL	5.8841 mL	11.7682 mL
10 mM	0.5884 mL	2.942 mL	5.8841 mL
50 mM	0.1177 mL	0.5884 mL	1.1768 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481