


4-Cyanobenzeneboronic acid

Chemical Properties

CAS No. :	126747-14-6
Formula:	C ₇ H ₆ BN ₂ O ₂
Molecular Weight:	146.939
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4-Cyanobenzeneboronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.8055 mL	34.0275 mL	68.055 mL
5 mM	1.3611 mL	6.8055 mL	13.611 mL
10 mM	0.6805 mL	3.4027 mL	6.8055 mL
50 mM	0.1361 mL	0.6805 mL	1.3611 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481