

Dipivaloylmethane

Chemical Properties


CAS No. : 1118-71-4

Formula: C₁₁H₂₀O₂

Molecular Weight: 184.2753

Storage: Keep away from direct sunlight, Store under nitrogen
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dipivaloylmethane is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4265 mL	27.1326 mL	54.2652 mL
5 mM	1.0853 mL	5.4265 mL	10.853 mL
10 mM	0.5427 mL	2.7133 mL	5.4265 mL
50 mM	0.1085 mL	0.5427 mL	1.0853 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481