


Boc-Propargylamine

Chemical Properties

| | |
|-------------------|--|
| CAS No. : | 92136-39-5 |
| Formula: | C ₈ H ₁₃ NO ₂ |
| Molecular Weight: | 155.1943 |
| Storage: | Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small> |



Biological Description

| | |
|-------------|--|
| Description | Boc-Propargylamine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 6.4437 mL | 32.2186 mL | 64.4371 mL |
| 5 mM | 1.2887 mL | 6.4437 mL | 12.8874 mL |
| 10 mM | 0.6444 mL | 3.2219 mL | 6.4437 mL |
| 50 mM | 0.1289 mL | 0.6444 mL | 1.2887 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481