


## Fmoc-OSu

## Chemical Properties

CAS No. :	82911-69-1
Formula:	C <sub>19</sub> H <sub>15</sub> N <sub>1</sub> O <sub>5</sub>
Molecular Weight:	337.3261
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Fmoc-OSu is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9645 mL	14.8223 mL	29.6446 mL
5 mM	0.5929 mL	2.9645 mL	5.9289 mL
10 mM	0.2964 mL	1.4822 mL	2.9645 mL
50 mM	0.0593 mL	0.2964 mL	0.5929 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481