


## 2-Bromo-5-iodopyridine

## Chemical Properties

CAS No. :	73290-22-9
Formula:	C <sub>5</sub> H <sub>3</sub> BrIN
Molecular Weight:	283.892
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	2-Bromo-5-iodopyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5225 mL	17.6125 mL	35.2249 mL
5 mM	0.7045 mL	3.5225 mL	7.045 mL
10 mM	0.3522 mL	1.7612 mL	3.5225 mL
50 mM	0.0704 mL	0.3522 mL	0.7045 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481