


2,4-Dichlorophenylboronic acid

Chemical Properties

CAS No. :	68716-47-2
Formula:	C ₆ H ₅ BCl ₂ O ₂
Molecular Weight:	190.82
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,4-Dichlorophenylboronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2405 mL	26.2027 mL	52.4054 mL
5 mM	1.0481 mL	5.2405 mL	10.4811 mL
10 mM	0.5241 mL	2.6203 mL	5.2405 mL
50 mM	0.1048 mL	0.5241 mL	1.0481 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481