


## (S)-N-Methylpyrrolidine-2-methanol

## Chemical Properties

CAS No. :	34381-71-0
Formula:	C <sub>6</sub> H <sub>13</sub> NO
Molecular Weight:	115.1735
Storage:	Pure form: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-N-Methylpyrrolidine-2-methanol is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.6828 mL	43.4141 mL	86.8282 mL
5 mM	1.7366 mL	8.6828 mL	17.3656 mL
10 mM	0.8683 mL	4.3414 mL	8.6828 mL
50 mM	0.1737 mL	0.8683 mL	1.7366 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481