


## 4-(Hydroxymethyl)benzeneboronic acid

## Chemical Properties

CAS No. :	59016-93-2
Formula:	C <sub>7</sub> H <sub>9</sub> BO <sub>3</sub>
Molecular Weight:	151.956
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	4-(Hydroxymethyl)benzeneboronic acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.5807 mL	32.9034 mL	65.8068 mL
5 mM	1.3161 mL	6.5807 mL	13.1614 mL
10 mM	0.6581 mL	3.2903 mL	6.5807 mL
50 mM	0.1316 mL	0.6581 mL	1.3161 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481