


## 2-Chloro-4-hydroxypyridine

## Chemical Properties

CAS No. :	17368-12-6
Formula:	C <sub>5</sub> H <sub>4</sub> ClNO
Molecular Weight:	129.544
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	2-Chloro-4-hydroxypyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.7196 mL	38.5981 mL	77.1962 mL
5 mM	1.5439 mL	7.7196 mL	15.4392 mL
10 mM	0.772 mL	3.8598 mL	7.7196 mL
50 mM	0.1544 mL	0.772 mL	1.5439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481