

Trans-2-octen-1-ol

Chemical Properties


CAS No. : 18409-17-1

Formula: C₈H₁₆O

Molecular Weight: 128.21

Storage: Store at RT

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Trans-2-octen-1-ol is an unsaturated alcohol compound widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.7997 mL	38.9985 mL	77.997 mL
5 mM	1.5599 mL	7.7997 mL	15.5994 mL
10 mM	0.780 mL	3.8999 mL	7.7997 mL
50 mM	0.156 mL	0.780 mL	1.5599 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wu L, et al. Geographical Indication Characteristics of Aroma and Phenolic Acids of the Changping Strawberry. Foods. 2023 Oct 24;12(21):3889.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481