

## trans,trans-2,4-Nonadienal

## Chemical Properties

CAS No. : 5910-87-2

Formula: C<sub>9</sub>H<sub>14</sub>O

Molecular Weight: 138.2069

Store under nitrogen

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	trans,trans-2,4-Nonadienal is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.2354 mL	36.1768 mL	72.3537 mL
5 mM	1.4471 mL	7.2354 mL	14.4707 mL
10 mM	0.7235 mL	3.6177 mL	7.2354 mL
50 mM	0.1447 mL	0.7235 mL	1.4471 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481