


2-Chloro-4-iodopyridine

Chemical Properties

CAS No. :	153034-86-7
Formula:	C ₅ H ₃ ClIN
Molecular Weight:	239.441
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-Chloro-4-iodopyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1764 mL	20.8821 mL	41.7641 mL
5 mM	0.8353 mL	4.1764 mL	8.3528 mL
10 mM	0.4176 mL	2.0882 mL	4.1764 mL
50 mM	0.0835 mL	0.4176 mL	0.8353 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481