


## 3-Ethynylpyridine

## Chemical Properties

CAS No. :	2510-23-8
Formula:	C <sub>7</sub> H <sub>5</sub> N
Molecular Weight:	103.1213
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	3-Ethynylpyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	9.6974 mL	48.4872 mL	96.9744 mL
5 mM	1.9395 mL	9.6974 mL	19.3949 mL
10 mM	0.9697 mL	4.8487 mL	9.6974 mL
50 mM	0.1939 mL	0.9697 mL	1.9395 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481