


Bis(benzonitrile)palladium chloride

Chemical Properties

CAS No. :	14220-64-5
Formula:	C ₁₄ H ₁₀ Cl ₂ N ₂ Pd
Molecular Weight:	383.57
Storage:	Keep away from direct sunlight, Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Bis(benzonitrile)palladium chloride is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6071 mL	13.0354 mL	26.0709 mL
5 mM	0.5214 mL	2.6071 mL	5.2142 mL
10 mM	0.2607 mL	1.3035 mL	2.6071 mL
50 mM	0.0521 mL	0.2607 mL	0.5214 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481