


## 2,3,5-Trichloropyridine

## Chemical Properties

CAS No. :	16063-70-0
Formula:	C <sub>5</sub> H <sub>2</sub> Cl <sub>3</sub> N
Molecular Weight:	182.435
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	2,3,5-Trichloropyridine is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4813 mL	27.4063 mL	54.8125 mL
5 mM	1.0963 mL	5.4813 mL	10.9625 mL
10 mM	0.5481 mL	2.7406 mL	5.4813 mL
50 mM	0.1096 mL	0.5481 mL	1.0963 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481