


4,6-Dichloropyrimidin-5-amine

Chemical Properties

CAS No. :	5413-85-4
Formula:	C ₄ H ₃ Cl ₂ N ₃
Molecular Weight:	163.993
Storage:	Keep away from moisture, Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4,6-Dichloropyrimidin-5-amine and dimethyl (5-pinol) isophthalate can be used to prepare an organic ligand through a Suzuki cross-coupling reaction.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0979 mL	30.4897 mL	60.9793 mL
5 mM	1.2196 mL	6.0979 mL	12.1959 mL
10 mM	0.6098 mL	3.049 mL	6.0979 mL
50 mM	0.122 mL	0.6098 mL	1.2196 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481