


(E)-Pent-2-enal

Chemical Properties

CAS No. :	1576-87-0
Formula:	C ₅ H ₈ O
Molecular Weight:	84.12
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	(E)-Pent-2-enal is an unsaturated aldehyde compound widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 80 mg/mL (951.02 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	11.8878 mL	59.4389 mL	118.8778 mL
5 mM	2.3776 mL	11.8878 mL	23.7756 mL
10 mM	1.1888 mL	5.9439 mL	11.8878 mL
50 mM	0.2378 mL	1.1888 mL	2.3776 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481