

## 2-Methyl-1-tetralone

## Chemical Properties


CAS No. : 1590-08-5

Formula: C<sub>11</sub>H<sub>12</sub>O

Molecular Weight: 160.2124

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Methyl-1-tetralone is a biochemical reagent utilized as a biomaterial or organic compound related to life sciences research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.2418 mL	31.209 mL	62.4181 mL
5 mM	1.2484 mL	6.2418 mL	12.4836 mL
10 mM	0.6242 mL	3.1209 mL	6.2418 mL
50 mM	0.1248 mL	0.6242 mL	1.2484 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481