


## Methyl 2,3-O-isopropylidene-5-O-tosyl-D-ribonucleoside

### Chemical Properties

CAS No. :	4137-56-8
Formula:	C <sub>16</sub> H <sub>22</sub> O <sub>7</sub> S
Molecular Weight:	358.407
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



### Biological Description

Description	Methyl 2,3-O-isopropylidene-5-O-tosyl-D-ribonucleoside is a biochemical reagent used in the synthesis of compounds.
-------------	---

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7901 mL	13.9505 mL	27.901 mL
5 mM	0.558 mL	2.7901 mL	5.5802 mL
10 mM	0.279 mL	1.3951 mL	2.7901 mL
50 mM	0.0558 mL	0.279 mL	0.558 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481