


## Butyramide

## Chemical Properties

CAS No. :	541-35-5
Formula:	C <sub>4</sub> H <sub>9</sub> NO
Molecular Weight:	87.12
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Butyramide is a fatty amide widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	11.4784 mL	57.3921 mL	114.7842 mL
5 mM	2.2957 mL	11.4784 mL	22.9568 mL
10 mM	1.1478 mL	5.7392 mL	11.4784 mL
50 mM	0.2296 mL	1.1478 mL	2.2957 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481