


## Fmoc-D-Trp-OH

## Chemical Properties

CAS No. : 86123-11-7  
 Formula: C<sub>26</sub>H<sub>22</sub>N<sub>2</sub>O<sub>4</sub>  
 Molecular Weight: 426.4639  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	Fmoc-D-Trp-OH (N $\alpha$ -Fmoc-D-tryptophan) is a biochemical reagent used as a biomaterial or organic compound in life science research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3449 mL	11.7244 mL	23.4489 mL
5 mM	0.469 mL	2.3449 mL	4.6898 mL
10 mM	0.2345 mL	1.1724 mL	2.3449 mL
50 mM	0.0469 mL	0.2345 mL	0.469 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481