


15-Crown-5

Chemical Properties

CAS No. :	33100-27-5
Formula:	C ₁₀ H ₂₀ O ₅
Molecular Weight:	220.2628
Storage:	Pure form: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	15-Crown-5 is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5401 mL	22.7004 mL	45.4009 mL
5 mM	0.908 mL	4.5401 mL	9.0802 mL
10 mM	0.454 mL	2.270 mL	4.5401 mL
50 mM	0.0908 mL	0.454 mL	0.908 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481