


Carboxybetaine methacrylate

Chemical Properties

CAS No. :	24249-95-4
Formula:	C11H19NO4
Molecular Weight:	229.2729
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Carboxybetaine methacrylate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3617 mL	21.8083 mL	43.6167 mL
5 mM	0.8723 mL	4.3617 mL	8.7233 mL
10 mM	0.4362 mL	2.1808 mL	4.3617 mL
50 mM	0.0872 mL	0.4362 mL	0.8723 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481