


2-Methylpropane-2-sulfinamide

Chemical Properties

CAS No. :	146374-27-8
Formula:	C ₄ H ₁₁ NOS
Molecular Weight:	121.201
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	2-Methylpropane-2-sulfinamide is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.2508 mL	41.2541 mL	82.5083 mL
5 mM	1.6502 mL	8.2508 mL	16.5017 mL
10 mM	0.8251 mL	4.1254 mL	8.2508 mL
50 mM	0.165 mL	0.8251 mL	1.6502 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481