


## 2-Bromo-3-pyridinol

## Chemical Properties

CAS No. :	6602-32-0
Formula:	C <sub>5</sub> H <sub>4</sub> BrNO
Molecular Weight:	173.995
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-Bromo-3-pyridinol is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7471 mL	28.7356 mL	57.4713 mL
5 mM	1.1494 mL	5.7471 mL	11.4943 mL
10 mM	0.5747 mL	2.8736 mL	5.7471 mL
50 mM	0.1149 mL	0.5747 mL	1.1494 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481