

## 5-Methoxyuridine 5'-triphosphate (trisodium) solution (100mM)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture,Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	5-Methoxyuridine 5'-triphosphate trisodium solution (100mM) is a modified nucleotide triphosphate (NTP). It can be incorporated into mRNA using T7 RNA polymerase.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481