


D-Glucurono-6,3-lactone acetonide

Chemical Properties

CAS No. :	20513-98-8
Formula:	C ₉ H ₁₂ O ₆
Molecular Weight:	216.188
Storage:	Store at RT <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	D-Glucurono-6,3-lactone acetonide can be transformed into optically active and partially protected inositol.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.6256 mL	23.1278 mL	46.2556 mL
5 mM	0.9251 mL	4.6256 mL	9.2511 mL
10 mM	0.4626 mL	2.3128 mL	4.6256 mL
50 mM	0.0925 mL	0.4626 mL	0.9251 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481