

5-Propargylamino-ddCTP (trisodium) solution (25mM)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 2 years, 4°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	5-Propargylamino-ddCTP (trisodium) solution (25mM) is a nucleoside molecule utilized in synthesizing anthocyanin dye-nucleotide conjugates, which are used for nucleic acid labeling or sequence analysis.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481