


Phenyl 2,3,4,6-tetra-O-acetyl- β -D-thioglucopyranoside

Chemical Properties

CAS No. :	23661-28-1
Formula:	C ₂₀ H ₂₄ O ₉ S
Molecular Weight:	440.464
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Phenyl 2,3,4,6-tetra-O-acetyl- β -D-thioglucopyranoside, 98% (PTATG, 98%), is suitable for glycoscience research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2704 mL	11.3518 mL	22.7035 mL
5 mM	0.4541 mL	2.2704 mL	4.5407 mL
10 mM	0.227 mL	1.1352 mL	2.2704 mL
50 mM	0.0454 mL	0.227 mL	0.4541 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481