


Acedoben sodium

Chemical Properties

CAS No. :	29305-16-6
Formula:	C ₉ H ₈ NNaO ₃
Molecular Weight:	201.1545
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Acedoben sodium is the sodium salt form of Acedoben. Acedoben can rapidly form a self-assembling coordination complex with iron ions, and this Fe-Ace coordination complex not only serves as a carrier for tumor antigens but also enhances antigen-specific antitumor immunity due to its intrinsic adjuvant properties.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9714 mL	24.8571 mL	49.7141 mL
5 mM	0.9943 mL	4.9714 mL	9.9428 mL
10 mM	0.4971 mL	2.4857 mL	4.9714 mL
50 mM	0.0994 mL	0.4971 mL	0.9943 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481