


## N-Cbz-1,2,3,6-tetrahydropyridine-4-boronic acid pinacol ester

## Chemical Properties

CAS No. :	286961-15-7
Formula:	C <sub>19</sub> H <sub>26</sub> BNO <sub>4</sub>
Molecular Weight:	343.225
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-Cbz-1,2,3,6-tetrahydropyridine-4-boronic acid pinacol ester is a biochemical reagent.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9135 mL	14.5675 mL	29.135 mL
5 mM	0.5827 mL	2.9135 mL	5.827 mL
10 mM	0.2913 mL	1.4567 mL	2.9135 mL
50 mM	0.0583 mL	0.2913 mL	0.5827 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481