


## Triamcinolone acetonide acetate

## Chemical Properties

CAS No. :	3870-07-3
Formula:	C <sub>26</sub> H <sub>33</sub> F <sub>0</sub> O <sub>7</sub>
Molecular Weight:	476.53
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Triamcinolone acetonide acetate (Triamcinolone acetonide 21-acetate) is a biochemical reagent utilized as a biomaterial or organic compound in life science research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0985 mL	10.4925 mL	20.985 mL
5 mM	0.4197 mL	2.0985 mL	4.197 mL
10 mM	0.2099 mL	1.0493 mL	2.0985 mL
50 mM	0.042 mL	0.2099 mL	0.4197 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481