


Lead carbonate(II)

Chemical Properties

CAS No. : 598-63-0
 Formula: CO₃Pb
 Molecular Weight: 267.2
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lead carbonate(II) is a biochemical reagent.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7425 mL	18.7126 mL	37.4251 mL
5 mM	0.7485 mL	3.7425 mL	7.485 mL
10 mM	0.3743 mL	1.8713 mL	3.7425 mL
50 mM	0.0749 mL	0.3743 mL	0.7485 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481