


## 5-Nitouracil

## Chemical Properties

CAS No. :	611-08-5
Formula:	C <sub>4</sub> H <sub>3</sub> N <sub>3</sub> O <sub>4</sub>
Molecular Weight:	157.0843
Storage:	Store at RT Actual storage temperature shall be subject to the COA.



## Biological Description

Description	5-Nitouracil is a biochemical reagent.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.3662 mL	31.8309 mL	63.6618 mL
5 mM	1.2732 mL	6.3662 mL	12.7324 mL
10 mM	0.6366 mL	3.1831 mL	6.3662 mL
50 mM	0.1273 mL	0.6366 mL	1.2732 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481