


Acetyl-DL-phenylglycine

Chemical Properties

CAS No. :	15962-46-6
Formula:	C ₁₀ H ₁₁ NO ₃
Molecular Weight:	193.1992
Storage:	Store at RT Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acetyl-DL-phenylglycine is a biochemical reagent.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.176 mL	25.8799 mL	51.7598 mL
5 mM	1.0352 mL	5.176 mL	10.352 mL
10 mM	0.5176 mL	2.588 mL	5.176 mL
50 mM	0.1035 mL	0.5176 mL	1.0352 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481