


15:0 Lyso PC

Chemical Properties

CAS No. :	108273-89-8
Formula:	C ₂₃ H ₄₈ N ₀ O ₇ P
Molecular Weight:	481.6035
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	15:0 Lyso PC is a biochemical reagent.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0764 mL	10.3821 mL	20.7641 mL
5 mM	0.4153 mL	2.0764 mL	4.1528 mL
10 mM	0.2076 mL	1.0382 mL	2.0764 mL
50 mM	0.0415 mL	0.2076 mL	0.4153 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481