


1,5-Anhydrosorbitol

Chemical Properties

CAS No. :	154-58-5
Formula:	C6H12O5
Molecular Weight:	164.16
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,5-Anhydrosorbitol is a validated marker of short-term glycemc control. This substance is derived mainly from food.
Targets(IC50)	Endogenous Metabolite

Solubility Information

Solubility	DMSO: 138.9 mg/mL (846.13 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Saline: 10 mg/mL (60.92 mM),Solution. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0916 mL	30.4581 mL	60.9162 mL
5 mM	1.2183 mL	6.0916 mL	12.1832 mL
10 mM	0.6092 mL	3.0458 mL	6.0916 mL
50 mM	0.1218 mL	0.6092 mL	1.2183 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Kathleen M Dungan. 1,5-anhydroglucitol (GlycoMark) as a Marker of Short-Term Glycemic Control and Glycemic Excursions. Expert Rev Mol Diagn. Jan 2008; 8 (1), 9-19

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481