


CCR2 antagonist 3

Chemical Properties

CAS No. :	1380100-86-6
Formula:	C ₁₇ H ₂₅ FN ₂ O ₂
Molecular Weight:	308.39
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	CCR2 antagonist 3 (AZD-2927) is an antagonist of CCR2.
Targets(IC50)	CCR

Solubility Information

Solubility	DMSO: 90 mg/mL (291.84 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 3.3 mg/mL (10.7 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2426 mL	16.2132 mL	32.4265 mL
5 mM	0.6485 mL	3.2426 mL	6.4853 mL
10 mM	0.3243 mL	1.6213 mL	3.2426 mL
50 mM	0.0649 mL	0.3243 mL	0.6485 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Jonas BOSTRÖM, et al. New ikach blockers. WO2012074469A1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481