


A33903

Chemical Properties

CAS No. : 70890-53-8
 Formula: C₁₇H₁₅N₃O₂
 Molecular Weight: 293.32
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	A33903 is a 1,4-benzodiazepine analog that is being studied as a potential inhibitor of respiratory syncytial virus. This molecule was moderately active and demonstrated good pharmacokinetic properties.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4092 mL	17.0462 mL	34.0925 mL
5 mM	0.6818 mL	3.4092 mL	6.8185 mL
10 mM	0.3409 mL	1.7046 mL	3.4092 mL
50 mM	0.0682 mL	0.3409 mL	0.6818 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481