

## Z-Hpg

## Chemical Properties


CAS No. : 1234692-80-8

Formula: C<sub>14</sub>H<sub>15</sub>NO<sub>4</sub>

Molecular Weight: 261.27

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Z-Hpg is a novel bortac-gj enabled selective probe for labeling of nascent proteomes, thus recording translation, in cell networks connected by gap junctions, leaving unconnected cells not labeled
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8275 mL	19.1373 mL	38.2746 mL
5 mM	0.7655 mL	3.8275 mL	7.6549 mL
10 mM	0.3827 mL	1.9137 mL	3.8275 mL
50 mM	0.0765 mL	0.3827 mL	0.7655 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481